


# Equine Gowns

## For Colic Surgeries

Our GEPCO Equine Colic Gowns come in a variety of styles and sizes, all of which feature heat-sealed seams that provide fluid protection between you and your patient. All GEPCO gowns are ready-to-use in convenient, sterile peel pouches to save you time, increase efficiency, and reduce potential for infection.




Size M - XL: PM.2094V

### **EQUINE COLIC GOWN - 2094V**

This gown is a great choice for those who desire good fluid protection with a gown that is cooler and more comfortable.

- ▶ Good fluid protection
- ▶ Strong, heavy-duty, yellow Tyvek® sleeves
- ▶ Sterile-back style
- ▶ Fluid-repellent, spunlace body
- ▶ Addition fluid-proof Tyvek® shield on inside front of gown
- ▶ Our most comfortable colic gown




Size M - XL: PM.2094VLSEP

### **EQUINE COLIC GOWN - 2094VLSEP**

This gown is a great choice for those who desire excellent fluid protection with a soft, smooth, and pliable sleeve.

- ▶ Excellent fluid protection
- ▶ Sturdy, blue, polyethylene sleeves
- ▶ Sterile-back style
- ▶ Impervious, laminated, spunlace body
- ▶ Slightly warm



Size M - XL: PM.2094VLS

### **EQUINE COLIC GOWN - 2094VLS**

This gown is a great choice for those who desire maximum fluid protection with a strong, rugged sleeve.

- ▶ Maximum fluid protection
- ▶ Strong, heavy-duty, yellow Tyvek® sleeves
- ▶ Sterile-back style
- ▶ Impervious, laminated, spunlace body
- ▶ Slightly warm

**Ask about our Customized Colic Surgery Packs!**